



# XELEX DNA SELEX CORE KIT



## PART 1: START OF SELECTION ROUNDS

KIT VERSION 1.0, MARCH 2013. SEE ALSO DETAILED VERSION OF THIS PROTOCOL.

A

### COUNTERSELECTION (OPTIONAL)

1



#### 1. Counterselection (optional step)

- Start with 1 - 3 µg purified dsDNA from the previous selection round. (First SELEX round: 10 - 50 µg ssDNA library).
- Mix dsDNA with 50 - 100 µl [10x] SELEX buffer (=1/10 final volume) and fill up with nuclease-free water to 500 - 1000 µl final volume.
- Denature dsDNA at 94°C for 3 min, then put immediately on ice.
- Add to magnetic beads without immobilized target structure.
- Incubate 1 h (round 1) or 15 min to 30 min (starting from round 2) at RT or 37°C.
- Capture magnetic beads with a magnet.
- Remove supernatant for use in the further selection process.

B

### APTAMER-TARGET BINDING

2



#### 2. Binding

- Immobilization: Bind target structure to beads (see separate protocol).
- Start with 1 - 3 µg purified dsDNA from counterselection step or from the previous selection round (First SELEX round: 10 - 50 µg ssDNA library).
- Mix DNA (in water or elution buffer) with 50 - 100 µl [10x] SELEX buffer (=1/10 final volume) and fill up with nuclease-free water to 500 - 1000 µl final volume.
- Denature DNA at 94°C for 3 min, then put immediately on ice (Round 2 and later).
- Add precooled DNA or ssDNA to magnetic beads with immobilized target structure
- Amount of beads: Round 1: 10fold- (max. 200 µl), later rounds: 1fold binding capacity amount.
- Add 500 µl 1x SELEX Buffer.
- Incubate approx. 1 h (round 1) or 15- 30 min (round 2 and later).
- Mix gently to prevent sedimentation of beads (gentle shaking or occasional pipetting).
- Incubation temperature: Room temperature for analyses, 37°C for in-vivo targets

C

### APTAMER-TARGET WASHING

3



#### 3. Washing

- Round 1: 1 x with 1 ml SELEX Buffer
- Round 2: 2 x with 500 µl - 1 ml SELEX Buffer
- Additional rounds: Variable, roughly one additional washing step per additional selection round.
- Varying and adjusting the selection stringency: Adjust DNA / RNA and magnetic beads amounts, binding time, number of washing steps, binding temperature and salt conditions.

D

### APTAMER-TARGET ELUTION

4



#### 4. Elution (final volume 100 µl, choose one of the following alternate methods)

- Direct PCR from DNA bound to beads
  - Introduce beads directly into PCR reactions without elution. Applicable for 1 µM beads only.
- Heat
  - Heat beads in dest. H<sub>2</sub>O + 2 mM EDTA to 70°C - 94°C for complete denaturation.
  - Fix beads to one side of tube using a magnetic device, quickly (!) pipet the supernatant off.
  - Transfer supernatant without any beads to a new reaction tube.
  - Precipitate DNA or purify on DNA spin columns.
- Heat + SDS (2 % [w/v])
  - Incubate beads in 2 % [w/v] SDS and heat to 70°C - 94°C.
  - Fix beads to one side of tube using a magnetic device, quickly (!) pipet the supernatant off.
  - Transfer supernatant without any beads to a new reaction tube.
  - Precipitate DNA or purify on DNA spin columns.
- Elute DNA with Binding Buffer.
  - Add 400 µl of orange-colored Orange-DX buffer and purify on DNA spin columns.
- Competitive elution with either natural ligand of immobilized target for specific displacement or excess amount of free target.
  - The latter approach may fail to enrich high affinity binders.
- For each selection round, keep backups of 50 % of (non-amplified) bead eluate.

THIS PROTOCOL IS AVAILABLE ONLINE

[WWW.ROBOKLON.COM](http://WWW.ROBOKLON.COM)



# XELEX DNA SELEX CORE KIT



## PART 2: EMULSION PCR PROTOCOL

KIT VERSION 1.0, MARCH 2013. SEE ALSO DETAILED VERSION OF THIS PROTOCOL.

### E

### EMULSION PCR - OIL PHASE PREPARATION



5

#### 5. Emulsion - Setting up the Oil Phase

- Prepare 300 µl Oil Surfactant Mixture per reaction (50 µl water phase). Mix  
220 µl Emulsion Component 1 (73 % [v/v])  
20 µl Emulsion Component 2 (7 % [v/v])  
60 µl Emulsion Component 3 (20 % [v/v])
- Mix thoroughly by vortexing.
- Keep on wet ice until further usage.

### F

### EMULSION PCR - WATER PHASE PREPARATION



6

#### 6. Nucleic Acid Amplification of Binders (choose one of the following alternate methods)

- - PCR / Taq DNA Polymerase (for DNA aptamer selection, protocol provided below)
  - Mutagenesis PCR (see separate protocol) (for DNA aptamer selection, ships with separate protocol which replaces step 6a)
  - NASBA / RNA Purification (for RNA aptamer selection, ships with separate protocol which replaces steps 6a - 14)

#### 6a. Emulsion - Setting up the PCR Water Phase

- SELEX PCR / Taq DNA Polymerase Protocol:  
Water Phase: Mix 50 µl PCR sample on ice.  
Template DNA: Introduce 50 - 130 pmol (= 1 - 3 µg) of purified dsDNA.
  - 10 x PCR Buffer (with or without MgCl<sub>2</sub>) 1 x 5 µl
  - MgCl<sub>2</sub> [25 mM], if added separately 1.5 mM (or 1-5 mM) 3 µl (or 2- 10 µl)
  - BSA, acetylated [10 mg/ml] 0.01 (or 0-1) mg/ml 0-5 µl
  - Bank40-5'-Primer [100 µM] 4 µM (or 0.2 - 2 µM) 2 µl
  - Bank40-3'-Primer [100 µM] 4 µM (or 0.2 - 2 µM) 2 µl
  - dNTP mix [5 mM] 400 µM 4 µl
  - Thermostable DNA Polymerase 2.5 U (or 1.25 - 2.5 U) 0.25 - 0.5 µl
  - Template DNA (bead eluate, or use beads with bound DNA as template) 1 - 3 µg, up to 10<sup>13</sup> copies 1-25 µl
  - Sterile, DNA-free H<sub>2</sub>O @50 µl

### G

### EMULSION PCR - MIXING OIL AND WATER PHASE



7

#### 7. Create Emulsion Reactions

- Mix 300 µl Oil Surfactant Mixture (precooled, 4°C) and 50 µl Water Phase.
- Vortex thoroughly for 5 minutes at 4°C on a vortexer with fixation aperture.
- Dispense emulsion to three thin-walled reaction tubes ("triplicates") and run ePCR.

8

#### 8. Emulsion PCR - Parameters

- Perform PCR reaction with the following parameters:  
95°C - 120 sec - 15 - 20 x (95°C - 30 sec, 55°C - 60 sec, 72°C - 180 sec) - 72°C - 300 sec

THIS PROTOCOL IS AVAILABLE ONLINE

[WWW.ROBOKLON.COM](http://WWW.ROBOKLON.COM)



# XELEX DNA SELEX CORE KIT



## PART 3: END OF SELECTION ROUNDS

KIT VERSION 1.0, MARCH 2013. SEE ALSO DETAILED VERSION OF THIS PROTOCOL.

### H

### DNA PURIFICATION - BINDING STEP

9

#### 9. Emulsion Breaking

- Apply 40 µl Activation Buffer SF to spin-column membrane, keep column at room temperature until usage, do not spin.
- Preheat Elution Buffer to 80°C (see section D, Elution).
- Pool triplicates of each sample into a 2 ml plastic reaction tube.
- Add 1 ml 2-butanol (or butanol) and break emulsion by vortexing.
- Add 500 µl of orange-colored buffer SFB (10 volumes of PCR water phase).
- Mix buffer completely with sample.
- Centrifuge for 2 min at 11 000 x g (approx. 12 000 rpm).

10

#### 10. DNA Binding Step

- Remove organic phase (leave a small remain, do not remove interphase).
- Transfer water phase, interphase and remains of organic phase into a spin-column / receiver tube assembly.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).

### I

### DNA PURIFICATION - WASHING STEPS

11

#### 11. First Washing Step

- Discard flow-through and place back spin-column.
- Add 600 µl of Wash-SFX buffer to spin-column.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).

12

#### 12. Second Washing Step

- Discard flow-through and place back spin-column.
- Add 350 µl of Wash-SFX buffer to spin-column.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).

13

#### 13. Removal of Wash Buffer Traces

- Discard flow-through and place back spin-column.
- Centrifuge for 2 min at 11 000 x g (approx. 12 000 rpm) to remove any remaining traces of Wash-SFX buffer.

### J

### DNA PURIFICATION - ELUTION STEP

14

#### 14. DNA Elution

- Place spin-column in new collection tube (1.5 - 2 ml).
- Add 50-150 µl Elution buffer (optional: heated to 80°C).
- Incubate for 2 min at room temperature.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).
- Discard spin-column, cap the collection tube.
- Determine DNA concentration and total PCR yield per reaction.
- If total PCR yield is lower or much lower than 0.5 µg:  
Use the largest possible volume of purified DNA as template for scale-up emulsion PCR (20 PCR cycles, restart from step 5).
- If total PCR yield equals or is higher than 0.5 (better more than 1 µg):  
Use 10-50% of DNA for new selection round and keep 50-90% of DNA as backup.
- Go to step 1 or 2 to start new selection round (total 3 - 15 rounds) or end SELEX and start library / aptamer analysis (steps 14 - 17).





# XELEX DNA SELEX CORE KIT



## PART 4: LIBRARY AND APTAMER ANALYSIS

KIT VERSION 1.0, MARCH 2013. SEE ALSO DETAILED VERSION OF THIS PROTOCOL.

### K

### DIVERSITY STANDARD PCR

15

#### 15. Diversity Standard PCR amplification (non-emulsified PCR)

- 10 x PCR Buffer B (with 1.5 mM MgCl<sub>2</sub>) 1 x 10 µl
- Bank40-5'-Primer [100 µM] 1 µM 1 µl
- Bank40-3'-Primer [100 µM] 1 µM 1 µl
- dNTP mix [5 mM] 400 µM 8 µl
- Thermostable DNA Polymerase [5U/µl] 2.5 U 0.5 µl
- Template DNA 0.5 µg, up to 10<sup>13</sup> copies @100 µl
- Sterile, DNA-free H<sub>2</sub>O



16

#### 16. Diversity Standard PCR - Parameters

- Perform non-emulsified PCR reaction with the following parameters:  
95°C - 120 sec - 2 - 4 x (95°C - 30 sec, 59°C - 30 sec, 72°C - 90 sec) - 72°C - 300 sec

### L

### DiVE ASSAY

17

#### 17a. DiVE Assay (alternate method to DiStRO assay)

- Transfer 200 ng aliquots of standard or sample to separate plastic reaction tubes. Duplicates: Four tubes per sample. Assay: S1 Nuclease treatment, control: untreated.
- Add [5x] S1-Nuclease buffer to assay and control tubes (final concentration: 1x). Mix thoroughly. Keep final volumes identical and as low as possible. Final volume : -----
- Denature and reanneal both the assay and the control tubes:
  - Denaturation: 3 min at 98°C
  - Reannealing: 5 min at 65°C
- Add 1 µl S1-Nuclease [1 U/µl] to the assay tube, but not to the control tube.
  - Incubate for 30 min at 65°C.
- Add 2 µl EDTA [0.5 M, pH 8.0] (final concentration: 2 mM).
- Analyse 20-40 ng equivalent assay and control tubes side-by-side on 2-3 % agarose gels.



### M

### DiStRO ASSAY

17

#### 17b. DiStRO Assay (alternate method to DiVE assay)

- Mix the following components:
  - 300 ng of each diversity-standard or DNA-sample
  - + 1 µl DA-buffer (10x)
  - + 1 µl SYBR-Green (1:1000 in nuclease free H<sub>2</sub>O)
- Fill up to 10 µl with nuclease free H<sub>2</sub>O

Real-Time-Cycler Program:

- Denaturation: 2 min at 95°C
- Annealing: 180 min at 76°C (or 2°C less than calculated T<sub>m</sub>)
  - Take one measurement per minute.
- Remelting: Cool to 20°C, then heat by incrementing + 0,5°C per 7 sec to 98°C.
  - Take one measurement per increment, i.e. per 7 sec.



### N

### FLAA ASSAY (LIBRARY AND CLONE ANALYSIS)

18

#### 18. FLAA Assay

- Binder Assay: Bind biotinylated target to streptavidin-coated microtiterplate wells. Add a 10-fold excess of target in relation to the denoted binding capacity. Incubate 1 h at room temperature.
- Negative control: Add 2 µl Biotin [5 mM]. Incubate for 10 min at room temperature.
- Wash plate twice with 250 µl SELEX buffer [1x].
- Dilute an appropriate amount of binder DNA in 50 µl SELEX Buffer [1x]. Amount : -----
- Heat denature the DNA solution (95°C, 2 min) and immediately cool down on ice.
- Transfer the precooled DNA solution to the target coated microtiter plate. Incubate over night at 4°C or for one to two hours at room temperature.
- Discard the DNA supernatant.
- Wash plate immediately before measurement with 100 µl binding buffer [1x].
- Add 50 µl Oligreen (or Picogreen) (1:500 in binding buffer).
- Measure twice after 9 min. (or after 4 min.): Ex. 485 nm, em. 527 nm.

