



RESTRICTION ENZYMES - DOUBLE DIGESTION

Enzyme	Dilution Buffer #	BSA [100 µg/ml]	Proto-type	Specificity	Reaction Buffers				Reaction temp. [°C]	Inactivation temp. [°C]	Recommended 10x Reaction Buffer
					Low	Medium	High	Acet			
AccI	1	+ BSA	AccI	GT [^] MKAC	100	75	<25	100	37	80	Acet
AcvI	1	+ BSA	PmaCI	CAC [^] GTG	75	100	75	100	37	80	Acet
AluI	2	+ BSA	AluI	AG [^] CT	75	75	50	100	37	65	Acet
ApaI	1	+ BSA	ApaI	GGGCC [^] C	25	50	0	100	25	65	Acet
AvaI	1	+ BSA	AvaI	C [^] YCGRG	25	100	50	25	37	80	Medium
Ball	1	+ BSA	Ball	TGG [^] CCA	<25	25	<25	100	37	65	Acet
BamHI	BamHI	+ BSA	BamHI	G [^] GATCC	<25*	100*	100	75*	37	NA	High
BanII	1	+ BSA	HgiIII	GRGCY [^] C	75	25	50	100	37	65	Acet
BglI	3	+ BSA	BglI	GCCN ₂ [^] NGGC	50	75	100	50	37	65	High
BglII	2	+ BSA	BglII	A [^] GATCT	<25	75	100	25	37	NA	High
BsiHKCI	1	+ BSA	AvaI	C [^] YCGRG	<25*	100	100	50*	65	NA	Medium
BssHII	1	+ BSA	BsePI	G [^] C CGCGC	100*	100	100	100*	50	80	High
BstXI	1	+ BSA	BstXI	CCAN ₅ [^] NTGG	<25	100	100	50	50	65	High
BsuTUI	1	+ BSA	Clal	AT [^] CGAT	0	75	50	100	37	65	Acet
CviJI	-	NO BSA	CviJI	RG [^] CY	NR	NR	NR	NR	37	65	CviJI
CviJI*	-	NO BSA	CviJI*	G [^] C	NR	NR	NR	NR	37	65	CviJI*
DpnI	-	NO BSA	DpnI	GA [^] TC	100	100	75	100	37	80	Acet
DraI	1	+ BSA	AhaIII	TTT [^] AAA	100	100	75	100	37	65	Acet
EcoRI	3	BSA + det.	EcoRI	G [^] AATTC	25*	100	100	50	37	65	High
EcoRV	1	+ BSA	EcoRV	GAT [^] ATC	<25*	75*	100	75*	37	80	High
FokI	-	+ BSA	FokI	GGATC (9/13) [^]	NR	100	NR	NR	37	65	Medium
HaeIII	1	+ BSA	HaeIII	GG [^] CC	50	100	50	75	37	80	Medium
HincII	2	+ BSA	HindII	GTY [^] RAC	50	100	100	100	37	65	Medium
HindIII	2	+ BSA	HindIII	A [^] AGCTT	50*	100	<25	75*	37	65	Medium
HinfI	1	+ BSA	HinfI	G [^] ANTC	75	100	75	75	37	80	Medium
HpaI	1	+ BSA	HpaI	GTT [^] AAC	25	50	25	100	37	NA	Acet
HpaII	1	+ BSA	HpaII	C [^] CGG	100	75	25	75	37	65	Low
KpnI	1	NO BSA	KpnI	GGTAC [^] C	100	50	<25	50	37	NA	Low
MboI	1	+ BSA	MboI	[^] GATC	75	100	100	100	37	65	Medium
MboII	1	+ BSA	MboII	GAAGA (8/7) [^]	100	50	25	100	37	65	Low
MluI	1	+ BSA	MluI	A [^] C GCGT	25	75	100	50	37	65	High
MmeI	-	NO BSA	MmeI	TCCRAC (20/18) [^] or (21/19) [^]	NR	NR	NR	NR	37	80	MmeI
MnlI	2	+ BSA	MnlI	CCTC (7/6) [^]	75	100	50	75	37	65	Medium
MspI	1	+ BSA	HpaII	C [^] CGG	100	100	25	75	37	65	Medium
NarI	1	+ BSA	NarI	GG [^] CGCC	NR	NR	NR	NR	37	65	NarI
NcoI	1	+ BSA	NcoI	C [^] CATGG	75	75	75	100	37	65	Acet
NdeI	1	+ BSA	NdeI	CA [^] TATG	50	100	75	100	37	65	Acet
NheI	-	+ BSA	NheI	G [^] CTAGC	50	0	0	50	37	65	NheI
NotI	3	+ BSA	NotI	GC [^] GGCCGC	0	75	100	25	37	65	High
NruI	1	+ BSA	NruI	TCG [^] CGA	0	25	100	25	37	65	High
PinAI	PinAI	+ BSA	AgeI	A [^] C CGGT	50	25	<25	100	37	65	PinAI
PstI	1	+ BSA	PstI	CTGCA [^] G	75	75	100	50	37	80	High
PvuI	2	+ BSA	PvuI	CGAT [^] CG	<25	75	100	50	37	80	High
PvuII	1	+ BSA	PvuII	CAG [^] CTG	75*	100	75	75	37	NA	Medium
RsaI	1	+ BSA	RsaI	GT [^] AC	100	75	50	75	37	65	Low
RsrII	3	+ BSA	RsrII	CG [^] GWCCG	50	75	<25	100	37	65	Acet
SacI	1	NO BSA	SacI	GAGCT [^] C	100	50	<25	100	37	65	Low
SacII	1	NO BSA	SacII	CCGC [^] GG	50	75	75	100	37	65	Acet
SalI	1	+ BSA	SalI	G [^] TCGAC	<25	<25	100	<25	37	65	High
Sau3AI	1	+ BSA	MboI	[^] GATC	75	50	<25	75	37	65	Sau3AI
Scal	1	+ BSA	Scal	AGT [^] ACT	<25*	50*	100	<25*	37	80	High
SfiI	3	+ BSA	SfiI	GGCCN ₄ [^] NGGCC	<25	100	25	100	50	NA	Medium
SinI	1	NO BSA	Avall	G [^] GWCC	75	50	25	100	37	NA	Acet
SmaI	1	+ BSA	SmaI	CCC [^] GGG	0	0	0	100	25	65	Acet
SpeI	-	+ BSA	SpeI	A [^] CTAGT	75	100	50	75	37	65	Medium
SphI	1	+ BSA	SphI	GCATG [^] C	75	100	50	100	37	65	Medium
SspI	1	BSA + det.	SspI	AAT [^] ATT	50*	100	50	50	37	65	Medium
StuI	1	+ BSA	StuI	AGG [^] CCT	75	100	50	75	37	65	Medium
TaqI	2	+ BSA	TaqI	T [^] CGA	25	50	100	50	65	NA	High
TaqII	Taq II Stor. Bf	NO BSA	TaqII	GACCGA (11/9) [^]	NR	NR	NR	NR	70	NA	TaqII
TspDTI	-	NO BSA	TspDTI	ATGAA (11/9) [^]	NR	NR	NR	NR	70	NA	TspDTI
TspGWI	-	NO BSA	TspGWI	ACGGA (11/9) [^]	NR	NR	NR	NR	70	NA	TspGWI
Tth111I	2	+ BSA	Tth111I	GACN [^] NNGTC	100	25	25	100	65	NA	Low
XbaI	1	+ BSA	XbaI	T [^] CTAGA	25*	100	100	75	37	65	Medium
XhoI	1	+ BSA	XhoI	C [^] TCGAG	75	100	100	75	37	65	Medium

NR - buffer is not recommended, to see recommended buffer check product specifications.

NA - thermal inactivation is not available

* - enzyme exhibits star activity under certain conditions

Version 3.0 - Nov 2010