

# GENEFLO 625 DNA Ladder, ROX Labeled

## *GENEFLO 625 DNA Ladder*

**Economical, fluorescent-labeled DNA Ladders for genotyping and gene sequencing; perfect for use with ABI PRISM™ genetic analyzers.**

**Cat. No.****Size****Description:**

E3125-01	100 loadings
E3125-02	800 loadings

**Storage Conditions:**

Store at +4°C in the DARK

- The uniform spacing of the GENEFL0 fragments ensures precise sizing throughout a 50-625 bp range.
- Compatible with all fluorescent based separation systems including the ABI PRISM.
- Consist of 24 evenly spaced bands ranging from 50-625 bp, plus 2 doublets (129 bp and 429 bp) for software interpretation and orientation.
- Fragments are spaced at intervals of 25 bp providing the perfect tool for genotyping analysis.
- Can be used in either gel or capillary electrophoresis.
- Composed of 26 fragments of the following sizes: 50, 75, 100, 125, 129, 150, 175, 200, 225, 250, 275, 300, 325, 350, 375, 400, 425, 429, 450, 475, 500, 525, 550, 575, 600 and 625.
- Makes genotyping affordable and efficient.

**Advantages:**

- Broader range, up to 625 bp.
- Evenly spaced, no gaps.
- Ensures accurate and consistent sizing of fragments.
- No split peaks.
- Doublets for easy software orientation.
- Labeled with ROX.